

## Supplementary

**Table S1: List of oligos used for Combined bisulfite restriction analyses**

431 forward	5' -GTTGTTGTGTGGTTGGGTGGGT-3'
431 reverse	5' -CCCAAAAAACTAACTAAACTATTCC-3'
127 forward	5' -GGTTTTGTTAGGGAGTAGTAG-3'
127 reverse	5' -CAAAACTACATACAAATAATTCC-3'
432 forward	5' -TTTATTGGTTAGTTAAGATGG-3'
432 reverse	5' -AAAATCAATACAAACCAACTACC-3'
411 forward	5' -GTATTTGTGTGGTATTGGAG-3'
411 reverse	5' -AACTAAAAAAAAACTCATATCCAC-3'