


## Supporting information



**S2 Fig. Protein pattern of tomato extracts.** (A) SDS-PAGE and (B) Western-Blot analysis of native tomato protein extracts exemplarily shown for different commercially available tomato cultivars. SDS-PAGE was performed under reducing conditions. Coomassie Brilliant Blue G250 was used for protein staining. For Western blot analysis specific polyclonal Sola l 4-antibody was used. The 18 kDa band, corresponding to the native Sola l 4, is marked with an arrow. M: PageRuler Prestained Protein Ladder.