			Long					Short			
	Spec	Sen	\mathbf{FPR}	Q	MCC	Spec	Sen	FPR	Q	MCC	
		AntiBP2									
Secreted control	59.6	73.0	0.51			60.5	71.0	0.47			
Bioactive	63.6	48.7	0.27			64.3	52.9	0.29			
All				60.9	0.22				62.0	0.24	
	CAMP										
Secreted control	61.3	62.0	0.39			59.2	70.7	0.49			
Bioactive	61.5	60.8	0.38			63.5	51.2	0.29			
All				61.4	0.23				60.9	0.22	
	PeptideRanker										
Secreted control	77.9	51.4	0.15			70.2	72.5	0.31			
Bioactive	63.7	85.4	0.49			71.6	69.2	0.27			
All				68.4	0.39				70.9	0.42	

Table S6. Independent test set with control peptide set from secreted proteins

Comparison of PeptideRanker (measured at a threshold of 0.5), CAMP and AntiBP2 tested on the independent test set, with control peptides randomly selected from secreted proteins. AntiBP2 did not return predictions for 223 of the long and 394 of the short peptides. CAMP did not return predictions for 11 of the long and 7 of the short peptides.